



Luminaire Based Accessoires

LCA8007/05 ActiLume Mntg Clip H513 25 W

Mounting clip for the ActiLume G2 Indus Sensor H513 when the sensor is to be mounted in the luminaire housing.

Mounting clip for the ActiLume G2 Indus Sensor H513 when the sensor is to be mounted in the luminaire housing.

Product data

Product Data		Numerator – packs per outer box	
Full product code	871829172792700	Numerator – packs per outer box	24
Order product name	LCA8007/05 ActiLume Mntg Clip H513 25 W	Material no. (12NC)	913700358803
EAN/UPC – product	8718291727934	Net weight (piece)	10.000 g
Order code	72792700		
Numerator – quantity per pack	1		