



Projection TV lamp replacement modules Christie

MODULE CHRISTIE 03-900471-01P

Using original Philips manufactured replacement modules, all with patented high quality Philips UHP lamp inside, is the best lamp replacement solution for your TV set. Philips replacement modules protect customers from any counterfeit modules. Philips replacement modules ensure the quality of projection experience in terms of safety, light output, reliability and efficiency. The portfolio of TV replacement modules covers the majority of projection TV brands.

Product data

Operating and Electrical		Material Nr. (12NC)	
Power Consumption	200 W	Numerator - Quantity Per Pack	1
Product Data		EAN/UPC - Product/Case	50046677421750
Order product name	MODULE CHRISTIE 03-900471-01P	Numerator - Packs per outer box	4
Full product name	MODULE CHRISTIE 03-900471-01P	EAN/UPC - Case	50046677421750
Order code	914499987595		