



Residential Collection

# Flowerbed

## Flowerbed post anthracite 1x11W 230V

Flowerbed - lantern post - antracit - 1 pc - 11 W

### Product data

General Information	
Number of light sources	1 pc
Lamp holder	2G7
Protection class IEC	Safety class I
Family name	Flowerbed
Product type	lantern post
Lamp type	PLS4

Operating and Electrical	
Maximum lamp power	11 W


  

Approval and Application	
Ingress protection code	IP44 [ Wire-protected, splash-proof]

Product Data	
Full product code	169169316
Order product name	Flowerbed post anthracite 1x11W 230V
EAN/UPC - Product	8717943747146
Order code	915002543701
Numerator - Quantity Per Pack	1
Numerator - Packs per outer box	3
Material Nr. (12NC)	915002543701
Net Weight (Piece)	2.500 kg

## Flowerbed

### Dimensional drawing



Residential Collection pedestal

