



InfraRed Healthcare Heat Incandescent

R95 IR 100W E27 230V Red 1CT/25

Philips infrared lamps for healthcare and bodycare applications are designed for treating deep-seated muscular ailments and sports injuries. These incandescent reflector lamps are an excellent solution to provide localized heat treatment to relieve muscular pain. They can also be used to treat rheumatic ailments. This form of heat therapy has also been shown to speed the healing of different kinds of injuries such as sports injuries and non-infected wounds, in many cases providing rapid and effective pain relief. The benefits of this form of heat therapy are based on locally enhanced blood circulation in the skin caused by vasodilatory response. This results in an increased transport rate of metabolites and other essential biochemical compounds. Benefits are also gained by deeper penetration of heat, which provides a gentle, pleasant warming effect.

Product data

General Information	
Cap-Base	E27 [E27]
Operating Position	UNIVERSAL [Any or Universal (U)]
Main Application	Infrared Health
Nominal Lifetime (Nom)	300 h

Operating and Electrical	
Power (Rated) (Nom)	100 W
Voltage (Nom)	230 V


Controls and Dimming	
Dimmable	Yes

Mechanical and Housing	
Bulb Finish	Red
Bulb Material	Soft Glass

Product Data	
Full product code	871150014559840
Order product name	R95 IR 100W E27 230V Red 1CT/25
EAN/UPC - Product	8711500145598
Order code	14559840
Numerator - Quantity Per Pack	1
Numerator - Packs per outer box	25
Material Nr. (12NC)	923244244208
Net Weight (Piece)	53.300 g

InfraRed Healthcare Heat Incandescent

Dimensional drawing



R95 IR 100W E27 230V Red

Product	D	C (max)
R95 IR 100W E27 230V Red 1CT/25	95 mm	130 mm

