



# PLUS T8 High Output

## F96T8/TL841/HO/PL ALTO TG

Philips PLUS T8 High Output 8-foot fluorescent lamps reduce the impact on the environment by offering long life and low mercury.

### Product data

| General Information         |   |
|-----------------------------|---|
| Cap-Base                    | R17D [ R17d (RDC)]                                      |
| Life 3-Hr Start [Hrs]       | 24000 h   |
| Features                    | ALTO® Plus  |
| Footnotes Fluorescent/CFL 1 | This Bulb Meets US Federal Minimum Efficiency Standard. |

| Light Technical                    |                             |
|------------------------------------|-----------------------------|
| Color Code                         | TL841 [ CCT of 4100K (841)] |
| Luminous Flux (Nom)                | 8200 lm                     |
| Color Designation                  | TL841                       |
| Design Mean Lumens                 | 7625 lm                     |
| Correlated Color Temperature (Nom) | 4100 K                      |
| Color Rendering Index (Nom)        | 85                          |

| Operating and Electrical |      |
|--------------------------|------|
| Power (Rated) (Nom)      | 86 W |

| Mechanical and Housing |            |
|------------------------|------------|
| Cap-Base Information   | Green Base |

|                       |    |
|-----------------------|----|
| Nominal Length (Inch) | 96 |
|-----------------------|----|

| Approval and Application   |               |
|----------------------------|---------------|
| Energy Saving Product      | Energy Saving |
| Picogram Per Lumen Hour    | 24 pg/lm.h    |
| Mercury (Hg) Content (Nom) | 4.4 mg        |


| Product Data                    |                           |
|---------------------------------|---------------------------|
| Order product name              | F96T8/TL841/HO/PL ALTO TG |
| EAN/UPC - Product               | 046677409128              |
| Order code                      | 927874084107              |
| Numerator - Quantity Per Pack   | 1                         |
| Numerator - Packs per outer box | 25                        |
| Material Nr. (12NC)             | 927874084107              |
| Net Weight (Piece)              | 0.001 kg                  |

# PLUS T8 High Output

## Dimensional drawing

### Product

F96T8/TL841/HO/PL ALTO TG



F96T8 86W/841 HO ALTO Plus ES

