



Medical Therapy Jaundice PL-L/PL-S

PL-L 18W/52/4P 1CT/25

By emitting light almost entirely within the 400 to 500 nm bandwidth these Medical Therapy Jaundice PL-L lamps have no radiation from the short wave UVB waveband. They are therefore ideal for treating new born babies suffering from hyperbilirubinemia (neonatal jaundice) and Crigler-Najjar Syndrome (CNS). Moreover, the bandwidth of these lamps peak at the most effective treatment wavelength of 450 nm. This highly efficacious phototherapy treatment has eliminated the need for blood transfusions in almost all jaundiced infants. In addition, with the compact format of this PL-L lamp equipment manufacturers have more design freedom in developing their solutions. By emitting a full spectrum of a high color temperature, the PL-L /953 lamps are ideal for treating SAD.

Produkt Daten

Allgemeine Eigenschaften		Spannung (Nom)	
Socket	2G11 [2G11]		57 V
Lebensdauer bis 50 % Ausfall vorgeheizt (Nom)	8000 h	Mechanische Kenndaten	
Lichttechnische Daten		Sockelinformation	4 Pin (4P)
Farbcode	52	Zulassungen und Anwendungseigenschaften	
Lichtstrom (Nom)	330 lm	Quecksilbergehalt (Nom)	4,4 mg
Restlichtstrom 2000 Std. (Nom)	93 %	Produktdaten	
Restlichtstrom 5000 Std. (Nom)	88 %	Gesamt-Produktcode	872790080517800
Farbkoordinate X (Nom)	155	Bestell-Produktname	PL-L 18W/52/4P 1CT/25
Farbkoordinate Y (Nom)	75	EAN/UPC - Produkt	8727900805178
Elektrische Kenndaten		Bestellcode	80517800
Power (Rated) (Nom)	17.5 W	Anzahl pro Verpackung	1
Lampenstrom (Nom)	0,375 A	Anzahl pro Umverpackung	25


Medical Therapy Jaundice PL-L/PL-S

Material-Nr. (12NC)	927904105206
Nettogewicht (Einzelteil)	60,900 g

Hinweise

- Es ist sehr unwahrscheinlich, dass der Glasbruch einer Lampe negative Auswirkungen auf Ihre Gesundheit hat. Wenn es zu einem Glasbruch kommt, lüften Sie den Raum 30 Minuten lang und entfernen Sie die Splitter möglichst mit Handschuhen. Legen Sie die Splitter in eine verschließbaren Plastikbeutel und geben Sie ihn beim Werkstoffhof zum Recycling ab. Benutzen Sie keinen Staubsauger.

Abmessungsskizzen



Product	D1 (max)	D (max)	A (max)	B (max)	C (max)
PL-L 18W/52/4P 1CT/25	18 mm	39 mm	188,2 mm	214 mm	220,6 mm

PL-L 18W/52/4P 219mm

