



MSR Gold™ FastFit

MSR Gold 2000 FastFit 1CT/8

All lamps burn out eventually, but when this happens with single ended MSR Gold™ FastFit, it is replaced in seconds – thanks to the specially designed lamp base and lamp holder. This lamp provides a high beam intensity of pure, white light for a truly illuminating performance, while the gold-plated caps provide superior heat protection and prevent premature failure. P3 technology allows use in any position and at higher temperatures, further extending lamp life and consistency of high-quality light output. Also, because the FastFit design is applied to Philips Halogen lamps, switching between lamp technologies can be done quickly and easily. The Philips MSR GoldTM 1200 FastFit can be operated in a lamp wattage range between 800W and 1400W.

Warnings and Safety

- A lamp breaking is extremely unlikely to have any impact on your health. If a lamp breaks, ventilate the room for 30 minutes and remove the parts, preferably with gloves. Put them in a sealed plastic bag and take it to your local waste facilities for recycling. Do not use a vacuum cleaner.

Product data


General Information		Luminous Efficacy (rated) (Nom)	
Cap-Base	PGJX50 [PGJX50]	Color rendering index (CRI)	91
Operating Position	UNIVERSAL [Any or Universal (U)]	Arc Length O (Nom)	8.0 mm
Life to 50% Failures (Nom)	750 h	Operating and Electrical	
Light Technical		Power Consumption	2,000 W
Color Code	- [Not Specified]	Lamp Current (Nom)	15.2 A
Luminous Flux (Min)	162,000 lm	Controls and Dimming	
Luminous Flux	180,000 lm	Dimmable	Yes
Chromaticity Coordinate X (Nom)	321		
Chromaticity Coordinate Y (Nom)	332		
Correlated Color Temperature (Nom)	6050 K		

MSR Gold™ FastFit

Product Data	
Order product name	MSR Gold 2000 FastFit 1CT/8
Full product name	MSR Gold 2000 FastFit 1CT/8
Full product code	872790091949300
Local code description	MSRG2000FF
Order code	928176905114
Material Nr. (12NC)	928176905114

Local order code	MSRG2000FF
Numerator - Quantity Per Pack	1
EAN/UPC - Product/Case	8727900919493
Numerator - Packs per outer box	8
EAN/UPC - Case	8727900919509

Dimensional drawing



Product	D (max)	O	L (min)	L (max)	L	C (max)	F
MSR Gold 2000 FastFit 1CT/8	34 mm	8.0 mm	66 mm	68 mm	67 mm	134 mm	50 mm

