

Premium LED bulbs Vintage Filament

LEDClassic50W ST64E27 2000K GOLD D 1BC/4

Featuring a classic heritage design, the new Philips Vintage LED lamp family is reminiscent of familiar incandescent bulbs. Available in an array of giant shapes, candles and other familiar form factors, they provide a seamless fit across your space. When used as a statement piece, they add a sense of nostalgia and create a highly desirable atmosphere for socialization and relaxing.

Product data

General Information	
Cap-Base	E27 [E27]
Nominal lifetime	15,000 hour(s)
Switching Cycle	20,000
Lighting Technology	LED
EU RoHS compliant	Yes
Light Technical	
Color Code	820 [CCT of 2000K]
Luminous Flux	630 lm
Color Designation	Flame
Correlated Color Temperature (Nom)	2000 K
Luminous Efficacy (rated) (Nom)	90.00 lm/W
Color Consistency	<6
Color rendering index (CRI)	80
LLMF At End Of Nominal Lifetime (Nom)	70 %
Operating and Electrical	
Line Frequency	50 to 60 Hz
Input Frequency	50 to 60 Hz


Power Consumption	7 W
Lamp Current (Nom)	40 mA
Wattage Equivalent	50 W
Starting Time (Nom)	0.5 s
Warm-up time to 60% light	0.5 s
Power Factor (Fraction)	0.7
Voltage (Nom)	220-240 V
Controls and Dimming	
Dimmable	Yes
Mechanical and Housing	
Bulb Shape	ST64
Approval and Application	
Energy Consumption kWh/1000 h	7 kWh
EyeComfort	Yes

Premium LED bulbs Vintage Filament

Product Data	
Order product name	LEDClassic50W ST64E27 2000K GOLD D 1BC/4
Full product name	LEDClassic50W ST64E27 2000K GOLD D 1BC/4
Full product code	871869668832800
Order code	929001228942


Material Nr. (12NC)	929001228942
Numerator - Quantity Per Pack	1
EAN/UPC - Product/Case	8718696688328
Numerator - Packs per outer box	4
EAN/UPC - Case	8718696688335

Dimensional drawing




Product	D	C
LEDClassic50W ST64E27 2000K GOLD D 1BC/4	64 mm	148 mm

Photometric data



Spectral Power Distribution Colour - LEDClassic50W ST64E27 2000K GOLD D 1BC/4



General uniform lighting - LEDClassic50W ST64E27 2000K GOLD D 1BC/4

