



# Premium LED bulbs Vintage Filament

## 5ST19/VIN/822/E26/CL-A/DIM BC 1PK

Featuring a classic heritage design, the new Philips Vintage LED lamp family is reminiscent of familiar incandescent bulbs. Available in an array of giant shapes, candles and other familiar form factors, they provide a seamless fit across your space. When used as a statement piece, they add a sense of nostalgia and create a highly desirable atmosphere for socialization and relaxing.

### Product data


<b>General Information</b>		<b>Power Consumption</b>	
Cap-Base	E26 [Single Contact Medium Screw]	5 W	
Nominal lifetime	15,000 hour(s)	<b>Lamp Current (Nom)</b>	
Switching Cycle	20,000	50 mA	
Lighting Technology	LED	<b>Wattage Equivalent</b>	
EU RoHS compliant	Yes	40 W	
<b>Light Technical</b>		<b>Starting Time (Nom)</b>	
Color Code	822 [CCT of 2200K]	0.5 s	
Luminous Flux	350 lm	<b>Warm-up time to 60% light</b>	
Color Designation	Flame	0.5 s	
Correlated Color Temperature (Nom)	2200 K	<b>Power Factor (Fraction)</b>	
Luminous Efficacy (rated) (Nom)	70.00 lm/W	0.7	
Color Consistency	<6	<b>Voltage (Nom)</b>	
Color rendering index (CRI)	80	120 V	
LLMF At End Of Nominal Lifetime (Nom)	70 %	<b>Temperature</b>	
<b>Operating and Electrical</b>		<b>T-Case Maximum (Nom)</b>	
Line Frequency	60 Hz	52 °C	
Input Frequency	60 Hz	<b>Controls and Dimming</b>	
		<b>Dimmable</b>	
		Yes	
		<b>Mechanical and Housing</b>	
		<b>Bulb Finish</b>	
		Amber	
		<b>Bulb Shape</b>	
		ST19	
		<b>Product Data</b>	
		<b>Order product name</b>	
		5ST19/VIN/822/E26/CL-A/DIM BC 1PK	

Premium LED bulbs Vintage Filament

Full product name	5ST19/VIN/822/E26/CL-A/DIM BC 1PK
Order code	929001334403
Material Nr. (12NC)	929001334403
Numerator - Quantity Per Pack	1
EAN/UPC - Product/Case	046677470401

Numerator - Packs per outer box	4
EAN/UPC - Case	50046677470406

Dimensional drawing



Product	D	C
5ST19/VIN/822/E26/CL-A/DIM BC 1PK	64 mm	138 mm

