



EvoKit Click 2x2

EvoKit CLKE 2x2 32L 24W 835 UNV SWZ P1

The Philips EvoKit LED Retrofit Kits are an energy efficient, easy to install solution to upgrade your fluorescent troffers to LED. Compatible with both standard and narrow T-grids, they offer a simple retrofit that will improve the look of your ceiling with its architectural styling without the need to actually break the ceiling plenum. The units also come standard with dimming capabilities, making them perfect for applications such as offices, classrooms, healthcare facilities, retail space and more. The Philips EvoKit LED Retrofit Kits offer the latest advances in LED technology, resulting in quality lighting with extremely high efficacies of up to 145 lm/w. The 1'x4', 2'x2' and the 2'x4' EvoKit are manufactured with quality components and finishes, meaning a consistent, balanced lighting scheme when using both configurations in the same space. The main diffuser and slanted troffer help reduce glare and create a pleasant, uniform throw of the light. Combine the aesthetics and quality with the ease of installation, and this product can literally transform your space in minutes!

Product data

General Information		Conduit Connection	
Lighting Technology	LED	7/8 inch knockout	
Mounting	Recessed	Mechanical and Housing	
Number Of Pieces	1	Shape	Square
Light Technical		Fixture Size	2 ft x 2 ft
Luminous Flux	3,200 lm	Enclosure	Metal and Plastic
Correlated Color Temperature (Nom)	3500 K	Product Data	
Luminous Efficacy (rated) (Nom)	133 lm/W	Order product name	EvoKit CLKE 2x2 32L 24W 835 UNV SWZ P1
Color rendering index (CRI)	80	Full product name	EvoKit CLKE 2x2 32L 24W 835 UNV SWZ P1
Operating and Electrical		Order code	EvoKit CLKE 2x2 32L 24W 835 UNV SWZ P1
Power Consumption	24 W	Material Nr. (12NC)	929001776313
Voltages	120/277 V	Numerator - Quantity Per Pack	1
		EAN/UPC - Product/Case	781087161906

EvoKit Click 2x2

Numerator - Packs per outer box	1
EAN/UPC - Case	50781087161901

