



MASTER LEDtube EM/ Mains T8

MASTER LEDtube Starter EMP GenII OP

The Philips MASTER LEDtube integrates a LED light source into a traditional fluorescent form factor. Its unique design creates a perfectly uniform visual appearance which cannot be distinguished from traditional fluorescent. For those that are looking for value for money within limited budget and re-lamping efforts for better light effect and lifetime.

Warnings and Safety

• -


Product data

General Information	
Lighting Technology	LED
CE mark	Yes
EU RoHS compliant	Yes
Product Data	
Order product name	MASTER LEDtube Starter EMP GenII OP
Full product name	MASTER LEDtube Starter EMP GenII OP
Full product code	871869681965400

Order code	81965400
Material Nr. (12NC)	929001801432
Numerator - Quantity Per Pack	1
EAN/UPC - Product/Case	8718696819654
Numerator - Packs per outer box	40
EAN/UPC - Case	8718696819678

MASTER LEDtube EM/Mains T8

Dimensional drawing



Product	D1	D2	C1
MASTER LEDtube Starter EMP GenII OP	3 mm	21.5 mm	34.5 mm

