



MASTER LEDtube EM/ Mains

MASTER LEDtube Starter EMP GenII OP

The Philips MASTER LEDtube integrates a LED light source into a traditional fluorescent form factor. Its unique design creates a perfectly uniform visual appearance which cannot be distinguished from traditional fluorescent. This product is the ideal solution for up lighting in general lighting applications.

Product data


Product Data		Numerator – packs per outer box	
Full product code	871869681965400		40
Order product name	MASTER LEDtube Starter EMP GenII OP	Material no. (12NC)	929001801432
EAN/UPC – product	8718696819654	Net weight (piece)	0.004 kg
Order code	81965400		
Numerator – quantity per pack	1		

Warnings and Safety

• -

MASTER LEDtube EM/Mains

Dimensional drawing



Product	D1	D2	C1
MASTER LEDtube Starter EMP GenII OP	3 mm	21.5 mm	34.5 mm

Dummy starter for cost down version

