



Lytecaster Accessory

Set of 4 Mounting Bar Clips

Set of 4 Mounting Bar Clips

Replacement cones for Lytecaster reflectors (not for use with Basic reflectors, Eyeballs and Shallow Specular reflectors, Pinholes and Slots). To order finishes other than Matte Black Step Baffle, add finish designations after catalog number (for example 1350CL for Specular Clear, 2000 Series).

Product data

General Information			
Mounting	Recessed	Material Nr. (12NC)	LIGU60000228
		Commercial Code	1956
		Catalog Number Description	Set of 4 Mounting Bar Clips
		Quick Ship	-
Product Data			
Order product name	Set of 4 Mounting Bar Clips		
Order code	1956		