



# Accessory Lens/Dichroic Color Filter

## 2" Dichroic Color Filter Medium Magenta

### 2" Dichroic Color Filter Medium Magenta

Color and diffusion filters, louvers and screens are available to precisely control and modify light for the most discerning spaces.

#### Product data

Product Data	
Order product name	2" Dichroic Color Filter Medium Magenta
Order code	ADF2MM
Material Nr. (12NC)	LIGU60002663
Commercial Code	ADF2MM

Catalog Number Description	2" Dichroic Color Filter Medium Magenta
Quick Ship	Spec Smart (10 Days) USA