



Accessory Lens/Dichroic Color Filter

2" Dichroic Color Filter Yellow Green

2" Dichroic Color Filter Yellow Green

Color and diffusion filters, louvers and screens are available to precisely control and modify light for the most discerning spaces.

Product data

Product Data	
Order product name	2" Dichroic Color Filter Yellow Green
Order code	ADF2YG
Material Nr. (12NC)	LIGU60002667
Commercial Code	ADF2YG

Catalog Number Description	2" Dichroic Color Filter Yellow Green
Quick Ship	-