



Calculite LED 6" cylinder gen 3

Calculite Stem Kit 36" White

Calculite Stem Kit 36" White

Calculite 6" Cylinder Gen3 features industry leading visual comfort and excellent uniform illumination over time.

Product data

General Information			
Lamp Source Code	LED	Commercial Code	CASK36WH
Downlight Type	Cylinders	Order code	CASK36WH
Mounting	Suspended	Catalog Number Description	Calculite Stem Kit 36" White
		Service / Spec Smart	-
		Material Nr. (12NC)	LIGU60014724
Product Data			
Order product name	Calculite Stem Kit 36" White		