



Renaissance LED pendant - large (RN20-RN30)

Renaissance, 80 LED, Type 4, Long Globe

Renaissance, 80 LED, Type 4, Long Globe - 116 lm/W

The Philips Lumec Renaissance LED pendant luminaires mix refinement together with ambition. The design reflects and evokes late 19th and early 20th century styling, perfectly suited for most urban and rural areas, while the state-of-the-art technology inside assures exceptional photometric performance, a long lifespan, and ease of maintenance. Designed to exceed your aesthetic expectations, the highly customizable Renaissance luminaires offer various decorative elements to add a feeling of sophistication to any area or project.

Product data

General Information		Initial Performance (IEC Compliant)	
Lamp Source Code	LED	Luminaire Lumens	15165
Color Temperature	4000K	Efficacy	116 lm/W
CRI	70	Product Data	
Photometric Distribution Type	Type IV	Order product name	Renaissance, 80 LED, Type 4, Long Globe
Mounting	Pole Arm	Commercial Code	RN-135W80LED4K-T-LE4R
Operating and Electrical		Order code	RN-135W80LED4K-T-LE4R
Luminaire Input Watts	130	Catalog Number Description	Renaissance, 80 LED, Type 4, Long Globe
Voltages	120/208/240/277/347/480 V	Service / Spec Smart	-
Mechanical and Housing		Material Nr. (12NC)	LUM70006882
Shape	Round		

Renaissance LED pendant - large (RN20-RN30)

