



Philips OneSpace luminous ceiling

Architectural, made-to-measure ceiling panel combining white LED light with textile

Architectural, made-to-measure ceiling panel combining white LED light with textile - 70 lm/W

Philips OneSpace luminous ceiling is an impressive custom sized (free hanging) unit with an excellent uniform light distribution, which contributes strongly to the acoustics in the room and adheres to the highest fire safety requirements. The result is a clutter free ceiling with an amazing light effect without seeing the light source.

Product data

General Information	
Color Temperature	4000K
CRI	80
Mounting	Suspended Surface Recessed
Operating and Electrical	
Voltages	120 V
Mechanical and Housing	
Shape	Rectangular
Initial Performance (IEC Compliant)	
Efficacy	70 lm/W

Product Data	
Order product name	Architectural, made-to-measure ceiling panel combining white LED light with textile
Commercial Code	OneSpace0001
Order code	OneSpace0001
Catalog Number Description	Architectural, made-to-measure ceiling panel combining white LED light with textile
Service / Spec Smart	-
Material Nr. (12NC)	PLC50000001